




## Complete Tamponade System

ebb is the only dual-balloon tamponade system designed to quickly and effectively treat postpartum hemorrhage.




See Package Insert for full instructions, warnings, contraindications, and precautions.

## ebb, Designed for Emergency Situations



### Conforms to Uterus

The larger (750mL) uterine balloon and malleable material allows ebb to conform to, and treat the uterus.



### Rapid Deployment

Simply spike an IV bag and gently squeeze to inflate the balloons in seconds, using the direct-fill ports.



### Safety Assurance



Monitor blood loss through the drain and verify lumen patency using the unique irrigation port.




### Eliminate Vaginal Packing

The dual-balloon design anchors the device, eliminating the need to pack the vagina and risk of retained pack.


~115,000 mothers with PPH in US each year<sup>5</sup>



### PPH RISK FACTORS<sup>3,6</sup>



### Clinical Evidence of ebb's Efficacy



**98%**  
ebb Significantly or Completely Stopped Bleeding

“the [ebb] is very useful in the management of PPH due to uterine atony and abnormal placentation”

**Avoided Hysterectomy 92%**

**45% >500mL**

45% patients required >500mL liquid to stop bleeding


**90%**

Vaginal Balloon Used

Data and quote:  
Initial experience with a dual-balloon catheter for the management of postpartum hemorrhage. Dildy GA, Belfort MA, Adair D, et al. Am J Obstet Gynecol. 2014; 210:136.e1-6.

### Product Support

ebb is complemented by an offering of training support for tamponade & postpartum hemorrhage.



#### Ordering Information

Product Description	Product Code	Box Quantity
ebb Complete Tamponade System	CTS-1000	1/ea
PPH Trainer (Training Accessory)	052-0083	1/ea

1. Clark SL, Belfort MA, Dildy GA. Maternal death in the 21st century: causes, prevention, and relationship to cesarean delivery. Am J Obstet Gynecol. 2008;199:36e1-36.e5.  
2. Pregnancy related mortality in the United States, 1998 to 2005. Berg CJ, Callaghan WM, Syverson C, Henderson Z. Obstet Gynecol. 2005;106:1228-1234.  
3. Prevention and management of postpartum haemorrhage. Royal College of Obstetricians and Gynaecologists (May 2009 with revisions April 2011). The Journal of Midwifery & Women's Health. 2004; 49(4):283-92.  
4. Obstetrics and Gynecology Chinese Medical Association. [Guideline of prevention and treatment about postpartum hemorrhage (2014)]. Zhonghua Fu Chan Ke Za Zhi 2014; 49: 641-646.  
5. Callaghan WM, Kuklina EV, Berg CJ. Trends in postpartum hemorrhage: United States, 1994 -2006. Am J Obstet Gynecol 2010;202:353.e1-6.  
6. Intrapartum care: care of healthy women and their babies during childbirth; NICE Clinical Guideline (Dec 2014).

Local Distributor:

**clinicalinnovations.com**  
747 West 4170 South | Murray, Utah USA 84123  
Phone 801.268.8200 | Fax 801.266.7373 |  
Toll Free 888.268.6222



**CLINICAL INNOVATIONS**  
forMOM. forBABY. forLIFE.